


Mark 95 & Mark 96 Bracket For Inlet Elbow Installation

The brackets enable users to mount a Mark 95 or Mark 96 regulator on a flat surface with enough elevation to allow the installation of a tri-clamp sanitary elbow (same size as inlet - no reducing elbows)



Part Number	Valve Size	Dimensions, Inches (mm)			
		A	B	C	D
95/96-150-LBKT	1-1/2"	8.83 (224)	2.00 (50,8)	5.40 (137)	0.63 (16,0)
95/96-2&3-LBKT	2"	11.07 (281)	2.00 (50,8)	6.75 (172)	0.63 (16,0)
95/96-2&3-LBKT	3"	11.07 (281)	2.00 (50,8)	7.75 (197)	0.63 (16,0)

How to Order:

Specify the part number above required for your particular application. Example: you need to mount a 2" MK95 on a flat structural support next to a WFI tank. In order to pipe the installation, the valve inlet must be connected to a 2" sanitary elbow. Select bracket model 95/96-2&3-LBKT.



Steriflow by Jordan Valve
 3170 Wasson Road • Cincinnati, OH 45209
 513.533.5600 • 800.543.7311 • 513.871.0105 (f)
 steriflow@richardsind.com • www.steriflowvalve.com